


# Actto 50



Included Components	QTY
Actto 50A Unit 220/230V	1
Power Cord	1
Reusable Electrode Set(6 type)	1
Dispersive Electrode with Cord	1
Footswitch assembly(Waterproof IP68)	1
Handpiece	1
Handpiece Stand	1
Owner's Guide	1

# Actto 50



## HIGH FREQUENCY THE COMPACT ELECTROSURGICAL UNIT

- + Exceptional cutting precision
- + Stop bleeding quickly
- + Fast healing

**DNH**medtech Co.

Address : #211/F-dong, 70 Jimok-ro, Paju-si, Gyeonggi-do, Korea(10880)

Tel : (+82)70-4242-1011 Fax : (+82)2-318-2727


e-Mail : DNHmedtech@naver.com URL : www.DNHmedtech.com

CE  
2460

DNH

**DNH MedTech provides technology with improved performance at low power settings minimizing the risk and damage, and stop bleeding quickly**

# Actto 50



## Actto 50 unit specifications

Width	10.4inches (265mm)
Height	3.5inches (89mm)
Depth	7.1inches (180mm)
Weight	6.8pounds (3.1kg)

### Output Power (10 Levels)

CUT	50Watts
CUT/COAG	40Watts
COAG	30Watts

Power Supply .... 220/230V, 50/60Hz, 220VA

Frequency ..... 1.2MHz ±10%

Electrical Classification ..... class1, BF type

Actto 50, a high frequency electrosurgery system helps manage soft tissues with greater efficiency, greater precision and less bleeding than the common dental scalpel

- Remove unwanted tissue with ease
- Reduce chair time by simplifying operative procedures and maintaining hemostasis
- Provides outstanding cutting precision for superior clinical and aesthetic results
- Allows uneventful healing

### Other clinical applications

- Widening gingival sulcus
- Crown lengthening
- Exposing impacted teeth
- Removing hyperplastic gingiva
- Frenectomy
- Stop bleeding quickly

## Clinical Cases


### • Gingivectomy

For extremely precise tracing of the ideal gingival line



### • Exposure of impacted tooth

Circular incision for a clean and successful exposure



### • Abscess incision and drainage

Fast and deep treatment



### • Dental neck exposure

Ideal for treating hidden decay



# Actto 50