

CS 3700

Scanning uncompromised.

Carestream | cs 3700

Design by STUDIO F·A·PORSCHE





High performance: A hallmark of Studio F. A. Porsche.





Revolutionary new workflow lets users navigate freely among different clinical applications, including restorative, orthodontic, implant and sleep device.

Gripping curves.

An ergonomically optimal grip channel empowers oral health professionals with a balanced sense of scanner control.

The added touch.

Touchscreen capability eliminates the extra steps of operation caused by switching back and forth between the scanner and the mouse.

Precision, accelerated.

The CS 3700 features turbo speed scanning that's tested and proven to scan in vitro single arches in 30 seconds, with the same proven accuracy as the CS 3600.

Strapping in.

The optional wrist strap connects user to the device for optimal control and care of the scanner.

Show your true colors.

Our patented Smart-shade color matching technology analyzes more data than other intraoral scanners, automatically detecting the enamel color of the scan area to identify the ideal match.

We approach color from all angles.

We don't determine the shade from a single RGB value. Instead, the CS 3700 features Bidirectional Reflectance Distribution Function (BRDF), which uses triangulation to factor in both the direction of illumination and the detection, providing more true-to-life color matching.



Shades of intelligence.

Unlike other systems, which collect color values from a single view, the CS 3700 collects actual shade values from several views of the tooth surface, leading to a more accurate match.

OTHERS: Texture-Based Shade Matching



CS 3700 BRDF



Smart-shade matching report.

In a single click, Smart-shade matching automatically generates a shade report for the selected area of interest, allowing users to easily share the outcome with the dental lab or referral.




Calibrated for brilliance.

10

The CS 3700 comes out of the box calibrated to provide the most reliable and stable shade matching results. Then, after 50 hours of usage, the smart-calibration tool alerts the user to calibrate the scanner through a simple 12-second guided calibration on a grey card.

Workflows that let you work freely.

CS ScanFlow software completely redefines intraoral scanning workflow. Unlike other scanning acquisition software, CS ScanFlow is completely open, letting you work freely among different clinical applications with unlimited workflow capacity, expanded processing options, one-click export and touchscreen capabilities.



Unlimited workflow capacity.

The completely new workflow design lets users navigate freely among different workflows from a single scan.

Expanded processing.

Expands data processing options, processing datasets quickly and accurately.

One-click export.

Simplifies data export and exchange.

Touchscreen capabilities.


An intuitive interface delivers direct software interaction via streamlined touchscreen navigation.

Hybrid scanning.

Enhances the capture of challenging data with hybrid scanning restorative and implant-borne restorative workflows.

Chairside Workflow

CS ScanFlow. Simplicity at work.



Shift your scanning into turbo.

The CS 3700 doesn't just provide fast scanning. It provides turbo speed scanning. In fact, this scanner has been tested and proven to scan in vitro single arches in 30 seconds when the recommended PC configuration is used. That's 20% faster than our previous scanner, without any compromise to the accuracy of the datasets.

TMO/ FASTER LU/O SCANNING

Processing that stays up to speed.

Turbo speed scanning requires turbo speed processing. The CS ScanFlow software lets you select from multiple processing options to ensure the best option for your indication, while giving you an average of 60% time savings on refinement.³



Backed by forward-thinking service.



The future is yours with **CS Advantage**. This premium service program provides a full array of service, support and warranty options for your purchase. For ultimate peace of mind, enroll on day one.



CS Protect

Extends the standard warranty and increases the length of coverage







© 2019 Carestream Dental LLC | 19094 AL 37 BR 0819 © 2019 Studio F•A•Porsche

1. When recommended PC configuration is used; in vitro test performed on a model.
2. Sending files to lab or third party via CS Connect requires enrollment in a services and support package.
3. With GPU (graphics processing unit) merging and improved texture mapping, the average time saving for refinement is approximately 60% on computers with 4 GB video memory.